# Professional Service Report の解説

## Professional Service Report の構成

TestingQA, SMTTest (A0000000)  
 Patient MRN: 000000 | DOB: JUN-01-1900 | Sex: Male  
 Diagnosis: Pancreatic ductal adenocarcinoma | Test Number 1

GUARDANT360<sup>CDx</sup>  
 Therapy Finder Page

REPORTING	PHYSICIAN
Report Date: JUN-01-20XX	XXXXXX XXXXX
Receipt Date: JUN-01-20XX	Account: XXXXXX Hospital
Collection Date: JUN-01-20XX	Address: XXXXXXXXXX, XXX-XXXX, Japan
Specimen: Blood	Ph: (XXX) XXXX-XXXX   Fax: (XXX) XXXX-XXXX
Status: FINAL	Additional Recipient: N/A

Complete Tumor Response Map on page 3

This content is provided as a professional service and has not been reviewed or approved by the FDA

### Summary of Detected Somatic Alterations, Immunotherapy Biomarkers & Associated Treatment Options

KEY: ✔ Approved in Indication ⚠ Approved in other indication ✘ Lack of response

Detected Alteration(s) / Biomarker(s)	Associated FDA-approved therapies	Clinical trial availability (see page 6)	% cfDNA or Amplification
MSI-High	<span style="color: green;">✔</span> Pembrolizumab <span style="color: orange;">⚠</span> Dostarlimab, Nivolumab, Nivolumab+ipilimumab	Yes	DETECTED
ATM S1455R	<span style="color: orange;">⚠</span> Olaparib, Talazoparib	Yes	24.2%
ATM K2811fs	<span style="color: orange;">⚠</span> Olaparib, Talazoparib	Yes	0.8%
BRCA2 T3033fs	<span style="color: orange;">⚠</span> Niraparib, Olaparib, Rucaparib, Talazoparib	Yes	0.7%
KRAS G12D	None	Yes	14.6%
GNAS R201H	None	Yes	10.1%
KRAS V14I	None	Yes	7.5%
CDKN2A R80*	None	Yes	6.1%
ARID1A R146I*	None	Yes	6.0%
TP53 R158H	None	Yes	0.1%
APC K2050fs	None	Yes	0.08%
BRAF Amplification	None	Yes	Medium (++)

This content is provided as a professional service and has not been reviewed or approved by the FDA

A more detailed Guardant360 Patient Report is available through our online portal: [portal.guardanthealth.com](http://portal.guardanthealth.com). To set up an account, contact Client Services: 855.698.8887

211-354-JP | RPT-002809 PB Pg 1 of 57

TestingQA, SMTTest (A0000000)  
 DOB: JUN-01-1900 | Test Number 1

GUARDANT360<sup>CDx</sup>  
 Tumor Biology Page

### Guardant360 Tumor Response Map

The Guardant360 laboratory developed test (LDT) Tumor Response Map illustrates the variant allele fraction (% cfDNA) of observed somatic variants at each sample submission. Amplifications are not plotted, and only the first and last five test dates are plotted. Please see the Physician Portal ([portal.guardanthealth.com](http://portal.guardanthealth.com)) for the Tumor Response Map with all test dates.

Highest Variant Allele Fraction: 24.6%

JUN 03, 2024

Detected Alteration(s) / Biomarker(s)	% cfDNA or Amp
ATM S1455R	24.2%
KRAS G12D	14.6%
GNAS R201H	10.1%
KRAS V14I	7.5%
HNF1A P291fs	7.5%
CDKN2A R80*	6.1%
ARID1A R146I*	6.0%
ATM K2811fs	0.8%
ERBB2 P1158P	0.6%
BRCA2 T3033fs	0.7%

This content is provided as a professional service and has not been reviewed or approved by the FDA

A more detailed Guardant360 Patient Report is available through our online portal: [portal.guardanthealth.com](http://portal.guardanthealth.com). To set up an account, contact Client Services: 855.698.8887

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**TestingQA, SMTTest (A0000000)**  
 Patient MRN: 000000 | DOB: JUN-01-1900 | Sex: Male  
 Diagnosis: Pancreatic ductal adenocarcinoma | Test Number 1

**GUARDANT360<sup>CDx</sup>**  
 Therapy Finder Page

---

**REPORTING**

Report Date: JUN-01-20XX  
 Receipt Date: JUN-01-20XX  
 Collection Date: JUN-01-20XX  
 Specimen: Blood  
 Status: FINAL

**PHYSICIAN**

XXXXX XXXXX  
 Account: XXXXX Hospital  
 Address: XXXXXXXXXXX, XXX-XXXX, Japan  
 Ph: (XXX) XXXX-XXXX | Fax: (XXX) XXXX-XXXX  
 Additional Recipient: N/A

Complete Tumor Response Map on page 3

---

This content is provided as a professional service and has not been reviewed or approved by the FDA

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Summary of Detected Somatic Alterations, Immunotherapy Biomarkers & Associated Treatment Options			
KEY <span style="color: green;">✔</span> Approved in indication <span style="color: orange;">~</span> Approved in other indication <span style="color: red;">✘</span> Lack of response			
Detected Alteration(s) / Biomarker(s)	Associated FDA-approved therapies	Clinical trial availability (see page 6)	% cfDNA or Amplification
MSI-High	<span style="color: green;">✔</span> Pembrolizumab <span style="color: orange;">~</span> Dostarlimab, Nivolumab, Nivolumab+ipilimumab	Yes	DETECTED
ATM S1455R	<span style="color: orange;">~</span> Olaparib, Talazoparib	Yes	24.6%
ATM K2811fs	<span style="color: orange;">~</span> Olaparib, Talazoparib	Yes	0.9%
BRCA2 T3033fs	<span style="color: orange;">~</span> Niraparib, Olaparib, Rucaparib, Talazoparib	Yes	0.6%
KRAS G12D	None	Yes	14.4%
GNAS R201H	None	Yes	10.5%

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**Variants of Uncertain Clinical Significance**  
 CDK6 R251M (6.3%), BRCA2 K385N (1.3%), APC S644G (0.9%), EGFR G503C (0.7%), CDK12 M840V (0.6%), CDK12 R890H (0.6%), RAF1 C95\* (0.4%), RET I657T (0.2%), CCND1 T12I (0.2%), ERBB2 A1232fs (Exon 27 deletion) (1.6%), ERBB2 L1131fs (Exon 26 deletion) (1.4%)  
 The functional consequences and/or clinical significance of alterations are unknown. Relevance of therapies targeting these alterations is uncertain.

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**Synonymous Alterations**  
 ERBB2 P1158P (0.7%), AR E566E (0.4%), TP53 S313S (0.4%), TERT R30R (0.3%), ERBB2 D1183D (0.2%)  
 This sequence change does not alter the amino acid at this position and is unlikely to be a therapeutic target. Clinical correlation is advised.

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Additional Biomarkers	
Biomarker	Additional Details
MSI-High	DETECTED

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
Detected Alteration(s) / Biomarker(s)	% ctDNA or Amp	
ATM S1455R	24.6%	
KRAS G12D	14.4%	
GNAS R201H	10.5%	
KRAS V14I	8.5%	
HNF1A P291fs	7.9%	
CDKN2A R80*	7.1%	
ARID1A R1461*	6.8%	
ATM K2811fs	0.9%	
ERBB2 P1158P	0.7%	Synonymous Alteration §
BRCA2 T3033fs	0.6%	
AR E566E	0.4%	Synonymous Alteration §
TP53 S313S	0.4%	Synonymous Alteration §
TERT R30R	0.3%	Synonymous Alteration §
TP53 R273C	0.3%	
ERBB2 D1183D	0.2%	Synonymous Alteration §
TP53 R158H	0.2%	
APC K2050fs	0.07%	
APC S644G	0.9%	Variants of Uncertain Clinical Significance §
BRAF Amplification Amplifications not graphed above	Medium (++) Plasma Copy Number 2.3	

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臨床試験 ID	試験名	実施地域		
<b>Available Clinical Trials (within the same state as the ordering physician)</b>				
There may be additional trials not listed here. Visit: <a href="http://portal.guardianhealth.com">portal.guardianhealth.com</a> or email <a href="mailto:clientservices@guardianhealth.com">clientservices@guardianhealth.com</a> with A0000000 in the subject line of the email, for additional trials.				
Alteration	Trial ID / Contact	Title	Phase	Site(s)
KRAS G12D	NCT05382559 Astellas Pharma Inc.,astellas. registration@astellas.com,800-888-7704	A Study of ASP3082 in Adults With Previously Treated Solid Tumors	Phase 1	Chuo-ku, Japan Koto-ku, Japan Kashiwa, Japan
	Visit <a href="http://portal.guardianhealth.com">portal.guardianhealth.com</a> for trials not within the same state as the physician's office			
CDKN2A R80*	NCT05999994 There may be multiple sites in this clinical trial. 1-877-CTLILLY (1-877-285-4559) or, Clinicaltrials.gov@lilly.com,1-317-615-4559	A Master Protocol (LY900023) That Includes Several Clinical Trials of Drugs for Children and Young Adults With Cancer	Phase 2	Chuo-ku, Japan Kobe, Japan Fukuoka, Japan
	Visit <a href="http://portal.guardianhealth.com">portal.guardianhealth.com</a> for trials not within the same state as the physician's office			
APC K2050fs	NCT04008797 Eisai Inquiry Service,eisai- chiken_hotline@hhc.eisai.co.jp	A Study of E7386 in Combination With Other Anticancer Drug in Participants With Solid Tumor	Phase 1	Nagoya, Japan Toyoake, Japan Hidaka, Japan Koto-ku, Japan  Additional trial sites available
	Visit <a href="http://portal.guardianhealth.com">portal.guardianhealth.com</a> for trials not within the same state as the physician's office			
MSI-High	Visit <a href="http://portal.guardianhealth.com">portal.guardianhealth.com</a> for trials not within the same state as the physician's office			
ATM S1455R	Visit <a href="http://portal.guardianhealth.com">portal.guardianhealth.com</a> for trials not within the same state as the physician's office			

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List of Available Clinical Trials				
Alteration	Trial ID / Contact	Title	Phase	Site (number in parenthesis is count of trial sites)
KRAS G12D	NCT04985604 Day One Biopharmaceuticals, clinicaltrials@dayonebio.com,650-484-0899	Tovorafenib (DAY101) Monotherapy or in Combination With Other Therapies for Patients With Melanoma and Other Solid Tumors	Phase 1 /Phase 2	Los Angeles, CA; Newport Beach, CA; Indianapolis, IN; Pittsburgh, PA; Portland, OR; Jacksonville, FL; Nashville, TN; Canada (2)
	NCT05111561 See <a href="https://clinicaltrials.gov/ct2/show/NCT05111561">https://clinicaltrials.gov/ct2/show/NCT05111561</a>	Testing the Combination of the Anticancer Drugs ZEN003694 and Binimetinib in Patients With Advanced/Metastatic or Unresectable Solid Tumors With RAS Alterations and Triple Negative Breast Cancer	Phase 1	Houston, TX; Boston, MA; Galveston, TX
	NCT05221320 Biomed Valley Discoveries, Inc., ERK@biomed-valley.com,816-960-6800	Trial of Ulixertinib in Combination With Hydroxychloroquine in Patients With Advanced Gastrointestinal (GI) Malignancies	Phase 2	Tucson, AZ; Saint Louis, MO; Cleveland, OH; Richmond, VA; New York, NY; New Brunswick, NJ; Fairway, KS; San Francisco, CA

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医薬品		薬剤の説明		現況
<b>Detailed Therapy Results</b>				
Alteration	Drug	Trade Name	Target	Current Status
MSI-High	ABL503		antiPD-L1/4-1BB bispecific antibody.	Phase 1 (Solid Tumor) Phase 1 (Brain and Central Nervous System Tumors)
	Adebrelimab		Anti-PD-L1 monoclonal antibody.	Phase 1 (Solid Tumor) Phase 3 (Small cell lung carcinoma (SCLC))
	AK129		Anti-PD-1/LAG-3 bispecific antibody.	Phase 1 (Solid Tumor) Phase 1 (Brain and Central Nervous System Tumors)
	AMG 404		Anti-PD-1 monoclonal antibody.	Phase 1 (Solid Tumor) Phase 1 (Brain and Central Nervous System Tumors)
	AP203		anti-CD137/PD-L1 bispecific antibody.	Phase 1 (Solid Tumor) Phase 2 (Non-small cell lung carcinoma (NSCLC), Head and neck squamous cell carcinoma (HNSCC), Esophageal squamous cell carcinoma)
	Atezolizumab	Tecentriq	Anti-PD-L1 monoclonal antibody.	Phase 2 (Pancreatic carcinoma) FDA Approved in other indications (Non-small cell lung carcinoma (NSCLC), Small cell lung carcinoma (SCLC), Alveolar soft part sarcoma)

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疾患における当該遺伝子の意義		薬剤感受性への影響	薬剤耐性との関連
<b>Relevance of Detected Alterations</b>			
Alteration	Role in Disease	Effect on Drug Sensitivity	Effect on Drug Resistance
<i>BRAF</i> V600E	<p>BRAF activating mutations or amplification have been reported to result in uncontrolled cell growth and tumorigenesis. <sup>(1,2)</sup> BRAF mutations have generally been reported to be mutually exclusive with RAS mutations in colorectal carcinoma (CRC), and have been associated with proximal tumor location, advanced tumor stage, poor tumor differentiation, and microsatellite instability. <sup>(3-7)</sup> A meta-analysis of 25 studies including 11,955 CRC patients reported that BRAF V600E mutations were associated with advanced TNM stage, poor differentiation, mucinous histology, microsatellite instability (MSI), CpG island methylator phenotype (CIMP), female gender, older age, proximal colon, and mutL homolog 1 (MLH1) methylation. <sup>(4)</sup></p>	<p>Braf signals upstream of the MAPK pathway, and BRAF amplification or activating mutations may confer sensitivity to inhibitors of Braf and/or components of the MAPK pathway, including MEK. <sup>(8)</sup> The BRAF V600-specific inhibitors vemurafenib and dabrafenib have been approved for the treatment of BRAF V600E-positive melanoma. <sup>(9,10)</sup> In addition, the MEK inhibitors trametinib and cobimetinib (in combination with vemurafenib) have been FDA-approved for BRAF V600E- and V600K-positive melanoma as has the encorafenib-binimetinib combination. <sup>(11-13)</sup> Vemurafenib has additionally been approved for BRAF V600-positive Erdheim-Chester disease. <sup>(14)</sup> Encorafenib in combination with cetuximab has been FDA-approved for the treatment of</p>	<p>In some studies, BRAF V600E mutation has been correlated with lack of response to approved anti-Egfr therapies such as cetuximab and panitumumab in CRC; patients with wild-type BRAF experienced a higher response rate when cetuximab was added to chemotherapy, compared to treatment with chemotherapy alone. <sup>(28,29)</sup> Several retrospective analyses (including the CRYSTAL and OPUS studies) suggest that while BRAF mutation is a negative prognostic marker, it is not predictive of response to cetuximab. <sup>(30-33)</sup></p>

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